



EMCO - DOPPLER

THE CASE

FOR

FETAL, ARTERIAL AND VASCULAR DIAGNOSIS



EMCO - DOPPLER MODEL D-520 For Easier Fetal Monitoring



A dedicated equipment designed for monitoring pregnancies using the principle of doppler ultrasound.

A complete self-contained system for fetal monitoring, with a permanently fixed 2 MHz transducer.

- Fetal heartbeats can be monitored from the 10th week of pregnancy
- Detection of multiple pregnancies.
- Localisation of placenta
- Identification of blood flow through umbilical cord
- Frequent monitoring throughout pregnancy provides constant reassurance to mother

Output can be presented in three different ways:

- Built-in speaker
- Stethoscope attachment

D520 (D)



Models Available ;
D520 (E) :- Economy (in Pouch)
D520(D) :- Deluxe (in Carrycase)

Compact, rugged, easy-to-use and reliable.

Hand held Battery powered audio unit

Auto power shut-off for maximum battery life.

Battery low indicator.

All accessories contained in a smart, tough attache-case. D520 D only

Can be connected to mains through A.C.Adaptor.