



### Pulse Oximeter Features

- 20 feet Fiber Optic Sensor (non-magnetic) avoids potentially hazardous heating or image artefacts during MR Scans
- Suitable for use in 0.5 to 3 Tesla MRI environment
- Built-In Rechargeable Battery back-up of 7 hours (Pulse Oximeter) & 6 hours (with NIBP)
- Alarm (Visual & Audio) for alarm limit violations

### NIBP Features

- Oscillometric method measuring Systolic, Diastolic, MAP & Pulse Rate
- Selectable Adult & Neonate functions with separate inflation and alarm settings
- Three modes MANUAL, STAT (5 quick measurements) & AUTO (measurements at user selectable time interval)
- Automatic pressure release for overpressure
- 20 ft Extension Tubing with Adult/Child/Infant Cuffs



 SunTech Medical<sup>®</sup>

Non-invasive Blood Pressure & Pulse Oximetry

Emco Meditek has been manufacturing MRI Compatible Pulse Oximeters since 2014, with an installation base all over India. Now comes the **MRI Compatible NIBP Oximeter 3700**, incorporating Nonin Medical (USA) SpO<sub>2</sub> & Suntech Medical (USA) NIBP Technology



At the heart  
of healthcare

**EMCO MEDITEK<sup>®</sup> PVT. LTD.**

An ISO9001 & ISO13485 Certified Company



